


Transurethral resection of the prostate/bladder (TUR-P/B) - monopolar



Materials/Instruments

- Resectoscope (optics, obturator, shaft, inner shaft, device insert with resection sling)
- Bladder syringe (glass bulb syringe)
- Twin tube for inflow to rinse the bladder
- Glycerol (ca. 20 litres)
- Compresses
- Mucous membrane desinfection and swabs
- Camera cover
- Electrical wire
- Adapter for tube
- Bin/bucket
- Gloves
- Irrigation catheter and syringe (10 ml)
- Lubricant (e.g. endosgel)

Electrical devices

- Electrical generator (Erbe, Autocon)
- Footswitch
- Monitor
- Camera
- Cold-light source

!!!After surgery it is important to maintain continuous bladder irrigation with sodium chloride!!!